



Sunrise Science Products

9925 Businesspark Ave., Suite B. San Diego, CA 92131

Phone: (760) 889-0982

Toll-Free Fax: (888) 880-9197

www.sunrisescience.com

PRODUCT DESCRIPTION

PRODUCT NAME:

SC Complete-CuSO₄ Powder

CATALOG NUMBER:

1489-XXX

STORAGE CONDITIONS:

Store tightly closed at 2-8°C

FINAL CONCENTRATION:

28.7 g/L

SUGGESTED USE:

Use at stated g/L. Filter sterilize or autoclave at 121°C for 15 minutes. To prepare plates, autoclave agar separately and add to sterile medium, or add agar to liquid medium and adjust to pH 5.8-6.0 before autoclaving.

CONTENTS/L:

YNB-CuSO₄ Powder (1.7 g/L)

Ammonium Sulfate (5 g/L)

Glucose (20 g/L)

Adenine (21 mg/L)

L-Alanine (85.6 mg/L)

L-Arginine (85.6 mg/L)

L-Asparagine (85.6 mg/L)

L-Aspartic Acid (85.6 mg/L)

L-Cysteine (85.6 mg/L)

L-Glutamine (85.6 mg/L)

L-Glutamic Acid (85.6 mg/L)

L-Glycine (85.6 mg/L)

L-Histidine (85.6 mg/L)

Inositol (85.6 mg/L)

L-Isoleucine (85.6 mg/L)

L-Leucine (173.4 mg/L)

L-Lysine (85.6 mg/L)

L-Methionine (85.6 mg/L)

Para-Amino Benzoic Acid (PABA) (8.6 mg/L)

L-Phenylalanine (85.6 mg/L)

L-Proline (85.6 mg/L)

L-Serine (85.5 mg/L)

L-Threonine (85.6 mg/L)

L-Tryptophan (85.6 mg/L)

L-Tyrosine (85.6 mg/L)

Uracil (85.6 mg/L)

L-Valine (85.6 mg/L)