



## **PRODUCT DESCRIPTION**

**PRODUCT NAME:**

SC-Glu-Gln Powder

**CATALOG NUMBER:**

6114-XXX

**STORAGE CONDITIONS:**

Store tightly closed at 2-8°C

**FINAL CONCENTRATION:**

1.83 g/L

**SUGGESTED USE:**

For complete media, add this amino acid mixture to YNB+Nitrogen+Sugar (dropout base) at stated g/L. Filter sterilize or autoclave at 121°C for 15 minutes. To prepare plates, autoclave agar separately and add to sterile medium, or add agar to liquid medium and adjust to pH 5.8-6.0 before autoclaving.

**CONTENTS/L:**

Adenine (21 mg/L)  
L-Alanine (85.6 mg/L)  
L-Arginine (85.6 mg/L)  
L-Asparagine (85.6 mg/L)  
L-Aspartic Acid (85.6 mg/L)  
L-Cysteine (85.6 mg/L)  
L-Glycine (85.6 mg/L)  
L-Histidine (85.6 mg/L)  
Inositol (85.6 mg/L)  
L-Isoleucine (85.6 mg/L)  
L-Leucine (173.4 mg/L)  
L-Lysine (85.6 mg/L)  
L-Methionine (85.6 mg/L)  
Para-Amino Benzoic Acid (PABA) (8.6 mg/L)  
L-Phenylalanine (85.6 mg/L)  
L-Proline (85.6 mg/L)  
L-Serine (85.5 mg/L)  
L-Threonine (85.6 mg/L)  
L-Tryptophan (85.6 mg/L)  
L-Tyrosine (85.6 mg/L)  
Uracil (85.6 mg/L)  
L-Valine (85.6 mg/L)