



## Sunrise Science Products

9925 Businesspark Ave., Suite B. San Diego, CA 92131

Phone: (760) 889-0982

Toll-Free Fax: (888) 880-9197

www.sunrisescience.com

### **PRODUCT DESCRIPTION**

#### **PRODUCT NAME:**

SC Complete-Leu Powder

#### **CATALOG NUMBER:**

1484-XXX

#### **STORAGE CONDITIONS:**

Store tightly closed at 2-8°C

#### **FINAL CONCENTRATION:**

28.53 g/L

#### **SUGGESTED USE:**

Use at stated g/L. Filter sterilize or autoclave at 121°C for 15 minutes. To prepare plates, autoclave agar separately and add to sterile medium, or add agar to liquid medium and adjust to pH 5.8-6.0 before autoclaving.

#### **CONTENTS/L:**

YNB Powder (1.7 g/L)  
Ammonium Sulfate (5 g/L)  
Glucose (20 g/L)  
Adenine (21 mg/L)  
L-Alanine (85.6 mg/L)  
L-Arginine (85.6 mg/L)  
L-Asparagine (85.6 mg/L)  
L-Aspartic Acid (85.6 mg/L)  
L-Cysteine (85.6 mg/L)  
L-Glutamine (85.6 mg/L)  
L-Glutamic Acid (85.6 mg/L)  
L-Glycine (85.6 mg/L)  
L-Histidine (85.6 mg/L)  
Inositol (85.6 mg/L)  
L-Isoleucine (85.6 mg/L)  
L-Lysine (85.6 mg/L)  
L-Methionine (85.6 mg/L)  
Para-Amino Benzoic Acid (PABA) (8.6 mg/L)  
L-Phenylalanine (85.6 mg/L)  
L-Proline (85.6 mg/L)  
L-Serine (85.5 mg/L)  
L-Threonine (85.6 mg/L)  
L-Tryptophan (85.6 mg/L)  
L-Tyrosine (85.6 mg/L)  
Uracil (85.6 mg/L)  
L-Valine (85.6 mg/L)