



Sunrise Science Products

127 Perimeter Park Rd., Suite A, Knoxville, TN 37922

Phone: (760) 889-0982

Toll-Free Fax: (888) 880-9197

www.sunrisescience.com

PRODUCT DESCRIPTION

PRODUCT NAME:

SC Complete-Trp-Ura Powder

CATALOG NUMBER:

1491-XXX

STORAGE CONDITIONS:

Store tightly closed at 2-8°C

FINAL CONCENTRATION:

28.62g/L

SUGGESTED USE:

Use at stated g/L. Filter sterilize or autoclave at 121°C for 15 minutes. To prepare plates, autoclave agar separately and add to sterile medium, or add agar to liquid medium and adjust to pH 5.8-6.0 before autoclaving.

CONTENTS/L:

YNB Powder (1.7 g/L)
Ammonium Sulfate (5 g/L)
Glucose (20 g/L)
Adenine (21 mg/L)
L-Alanine (85.6 mg/L)
L-Arginine (85.6 mg/L)
L-Asparagine (85.6 mg/L)
L-Aspartic Acid (85.6 mg/L)
L-Cysteine (85.6 mg/L)
L-Glutamine (85.6 mg/L)
L-Glutamic Acid (85.6 mg/L)
L-Glycine (85.6 mg/L)
L-Histidine (85.6 mg/L)
Inositol (85.6 mg/L)
L-Isoleucine (85.6 mg/L)
L-Leucine (173.4 mg/L)
L-Lysine (85.6 mg/L)
L-Methionine (85.6 mg/L)
Para-Amino Benzoic Acid (PABA) (8.6 mg/L)
L-Phenylalanine (85.6 mg/L)
L-Proline (85.6 mg/L)
L-Serine (85.5 mg/L)
L-Threonine (85.6 mg/L)
L-Tyrosine (85.6 mg/L)
L-Valine (85.6 mg/L)